




dimensions of the PTF MOSAIC Dewar flask:
 12"x12"x6.5"
 XY dimensions of the removable panel = 9.75"x4.25"

Cryotiger standard cold end (http://www.ociw.edu/instrumentation/ccd/imacs/docs/ICG_Brochure.pdf):



High performance cold end (http://www.ociw.edu/instrumentation/ccd/imacs/docs/ICG_Brochure.pdf):



Inc.

CCD package (Data sheet is located at: <http://www.ociw.edu/instrumentation/ccd/parts/CCID-20.pdf>)

